



**DECTRIS**

Next Generation X-Ray Detectors

# News about DECTRIS and PILATUS

Christian Broennimann

CEO

DECTRIS Ltd



**DECTRIS**

Next Generation X-Ray Detectors

## Outline

- ▶ About DECTRIS
- ▶ PILATUS 6M: Status and Results
- ▶ PILATUS 2M: Status and Results



**DECTRIS**

Next Generation X-Ray Detectors

## About DECTRIS

- ▶ Spin-off Company of the Paul Scherrer Institut (PSI)
- ▶ Produces and sells next generation of X-ray detectors
- ▶ Founded 28<sup>th</sup> September 2006 as a private limited company
- ▶ DECTRIS has an exclusive license
- ▶ DECTRIS remains closely connected to PSI to stay at the cutting edge of detector technology

## DECTRIS

- ▶ Christian Broennimann, Ph. D., CEO
- ▶ Eric F. Eikenberry, Ph. D., Software Development
- ▶ Markus Näf, Design and Production
- ▶ Petr Salficky, Marketing and Sales
  
- ▶ Miro Kobas, Software Development, Integration
- ▶ Taniel Sakhelashvili, Production
  
- ▶ 2-3 open positions



## Our Products




**PILATUS 100K**  
1 module  
Area: 83.8 x 33.5 mm<sup>2</sup>  
Frame Rate: 300Hz



**PILATUS 500K**  
1 x 5 modules  
Area: 33.5 x 431 mm<sup>2</sup>  
Frame Rate: 12Hz



**PILATUS 6M**  
5 x 12 modules  
Area: 431 x 448 mm<sup>2</sup>  
Frame Rate: 12Hz



**PILATUS 2M**  
3 x 8 modules  
Area: 254 x 289 mm<sup>2</sup>  
Frame Rate: 30Hz

## Our customers

### ▶ Synchrotron sources

- ▶ SLS (CH), SPring8 (Jpn), Diamond (UK), APS (US), Hasylab and EMBL Hamburg (Ge)

- ▶ 15 PILATUS 100K sold
- ▶ 7 PILATUS 100K delivered
- ▶ 1 PILATUS 500K sold and ready to be delivered
- ▶ 2 PILATUS 2M sold to Diamond. Delivery in 2. Quarter 2008

### ▶ Alcatraz C-Mod, MIT Boston (Fusion reactor)


- ▶ 4 PILATUS 100K

### ▶ Industrial XRD

- ▶ 5 PILATUS 100K

## PILATUS 6M: 9 years of development

- ▶ Based on hybrid pixel technology developed by PSI for high energy physics
- ▶ Technology adapted to meet requirements of synchrotron radiation research
- ▶ Developed specifically for macromolecular crystallography



## PILATUS 6M: Specifications

- ▶ No of Modules: 60
- ▶ Module arrangement: 5 x 12
- ▶ Detector Size: 431 x 448 mm<sup>2</sup>
- ▶ Format: 2527 x 2463 pixels
- ▶ Spatial resolution: 0.172 x 0.172 mm<sup>2</sup>
- ▶ Dynamic range/pixel: 20bits
- ▶ Count rate/pixel: <8 Mcps/pixel
- ▶ Readout time: 3.5 ms
- ▶ Frame rate: 12.5 Hz
- ▶ Mechanical positioning: Flat geometry
- ▶ Intermodule gap: x: 7 pixels, y: 17 pixels,  
8.4% of total area








**DECTRIS**

Next Generation X-Ray Detectors

## Specifications cont.

- ▶ Pixel yield: >99.99% ( $1.6 \times 10^{-4}$  unreliable pixels)
- ▶ Module alignment: < 1.5 pixel
- ▶ Lower level threshold: 4 -12 keV
- ▶ Threshold Dispersion: 50 eV
- ▶ Energy resolution: 500 – 600 eV
- ▶ Shutter synchronization via external trigger input
- ▶ Operation modes:
  - 1) Continuous mode, exposure time defined by detector, shutterless operation
  - 2) Shutter operated mode, exposure time externally defined
- ▶ Data formats: Raw data, TIFF, EDF, CBF (compressed)

## PILATUS 6M at X06SA






# PILATUS 6M Integration at X06SA@SLS

May 2007

- 4 GB Fibre Channel
- 1 GB Ethernet



**DAQ Server:**

- 8 Cores
- 32 GB RAM

**DAQ RAID:**


- 10 TB Storage
- RAID 5

**Computing Cluster**

- 8 x 4 = 32 Cores
- 8 x 8 GB = 64 GB RAM


**GPFS Fileserver**

- 10 TB per Storage
- RAID 6
- 40 TB total GPFS Storage



green:  $f_o - f_c$  density of A-site  
Blue:  $2f_o - f_c$  density of rRNA


30S ribosomal subunit, Ramakrishnan Lab



Blue:  $2f_o - f_c$  density of helix 44

## Cubic Insulin: Fine $\phi$ -slicing

angular speed = constant



-> Better R-factors with fine- $\phi$  slicing

XDS artifacts at very high resolution





## R-Factors

Resolution [Å]	1 ° 1 s	0.5 ° 0.5 s	0.25 ° 0.25 s	0.125 ° 0.125 s	0.1 ° 0.1 s
5	2.20%	2.10%	2.30%	1.80%	2.00%
4	2.00%	2.00%	1.70%	1.60%	2.00%
3	2.10%	1.90%	1.80%	1.80%	2.10%
2.5	2.60%	2.20%	2.10%	2.20%	2.30%
2.2	3.30%	2.70%	2.40%	2.50%	2.70%
2	4.00%	3.50%	3.20%	3.10%	3.40%
1.8	5.00%	4.30%	4.00%	3.90%	5.00%
1.5	9.20%	8.00%	7.50%	7.60%	10.10%
1.4	19.60%	16.70%	15.60%	15.90%	20.00%
1.3	36.90%	31.80%	25.10%	24.50%	28.20%
1.2	49.40%	42.80%	34.50%	29.10%	29.60%



## Software Support (October 2007)

- ▶ XDS, W. Kabsch, 3D profile fitting
- ▶ MOSFLM, A. Leslie, 2D profile fitting
- ▶ D\*Trek (J. Pflugrath, Rigaku)



Experimental maps  
based on D\*Trek  
and a SAD data-set



**DECTRIS**

Next Generation X-Ray Detectors

## Status at BL X06SA in September 2007

- ▶ Hardware and Software-Integration finished
- ▶ Regular user operation started in June 2007
- ▶ 35 user groups so far
- ▶ Processing with newest XDS version (W. Kabsch) and MOSFLM
- ▶ Usually similar or lower R-factors with fine- $\phi$  slicing achieved than with CCD, analysis still ongoing
- ▶ Detector system 4-5x faster than CCD




## User Feedback

- ▶ "Faster data collection"
- ▶ "Somewhat less radiation damage"
- ▶ all but 2 user groups report at least as good data as with CCD
- ▶ "We are the limiting factor"
- ▶ "We don't have time to smoke a cigarette or to drink coffee"

- Icosahedral AlCuFe quasicrystal
- energy threshold set to 8 keV or 10 keV

Fluorescence suppression improves data quality in

- ⇒ phasing experiments
- ⇒ diffuse scattering



## Pilatus 2M Specs

- ▶ No of Modules: 24
- ▶ Detector Size 254 x 289 mm<sup>2</sup>
- ▶ Number of modules 3 x 8 = 24
- ▶ Format 1475 x 1679 = 2'476'525 pixels
- ▶ Spatial resolution: 0.172 x 0.172 mm<sup>2</sup>
- ▶ Dynamic range/pixel: 20bits
- ▶ Count rate/pixel: ~3 MHz/pixel
- ▶ Readout time: 3.5 ms
- ▶ Frame rate: 30 Hz
- ▶ Power consumption 200 W
- ▶ Operating temperature 25°C
- ▶ Dimensions 388 x 434 x 526 mm<sup>3</sup>
- ▶ Weight Approx. 50 kg








**DECTRIS**

Next Generation X-Ray Detectors

## Pilatus 2M at BL X12SA




## Results : Ag behenate



Silver behenate, 12 keV beam, Threshold 55 %. Vrf -0.3 V. Sample-Detector ca. 2,2 m.

## Results : Collagen Sample



One out of 45000 images collected of a collagen bone probe



**DECTRIS**

Next Generation X-Ray Detectors

## Conclusions

- ▶ PILATUS 2M and 6M integration finished at SLS
- ▶ Break-through in speed and data-quality
- ▶ Data routinely collected in user operation mode
- ▶ Huge amount of data collected, storage and data processing issues should be emphasized
- ▶ DECTRIS is on schedule

**Thank you for your attention**